

UHP+

Ultrasound Hygiene Protocol Plus™

Your smile. Your comfort.
Your health.



Soniprox⁺

www.soniprox.com

UHP+ Daily hygiene care at its finest.



Only certified clinics guarantee the authenticity of the UHP+ protocol and your complete safety.

Soniprox⁺

www.soniprox.com

Your Daily UHP+™ Plan

Morning:

1. Place your aligners or denture in the Soniprox+ Ultrasound Cleaner with water and a cleaning tablet.
2. Select the 5-minute cleaning cycle.
3. Remove the appliance and rinse it under running water.
4. Put the attachments or denture back into your mouth.
5. Empty and rinse the tank, then set the UV cleaning program (Program 3) for 15 minutes. After completion, the device will turn off automatically.

During the day (up to 4×):

Use Soniprox+ CleaningFoam in two ways:

1. Method 1: Apply the foam directly onto the aligners and put them back in (spit out excess).
2. Method 2: Apply foam into your mouth, spread it with your tongue, spit out excess, then insert aligners.

Additionally, you can use the soft toothbrush from the Soniprox+ TravelSet to clean stubborn deposits on the delicate surfaces of the appliance. You can also take extra care of tongue hygiene by using the dedicated tongue-cleaning spatula included in the Soniprox+ TravelSet.

Evening:

1. Clean again with the 10-minute cycle in DEGAS mode.
2. Brush your teeth thoroughly and use interdental brushes for the best results.
3. Insert the aligners overnight. Store cleaned dentures overnight according to your doctor's recommendations.
4. Empty and rinse the tank of the Soniprox+ Ultrasound Cleaner, then set the UV cleaning program (Program 3) for 15 minutes. After completion, the device will turn off automatically.

Key rules:

- ✗ Do not use toothpaste on aligners.
- 🚫 Use only cold or lukewarm water.
- 🚫 Regularly disinfect the container and device.
- ✓ Use Soniprox+ products daily as part of therapy.
- ✓ Change the water in the device after each cleaning of the instruments.

Why is it worth it?

1. Protection against tooth decay and gum disease.
2. Prevention of inflammation and fungal infections.
3. Fresh breath and daily comfort.
4. Cleanliness and safety of orthodontic aligners, retainers, and dentures.
5. Support for orthodontic and prosthetic treatment.

Ekosystem SONIPROX+



1 Soniprox+ Ultrasonic Cleaner

- ✓ Provides deep, safe cleaning that protects against micro-damage.
- ✓ Removes biofilm, bacteria, stubborn deposits, and neutralizes odors.
- ✓ Extends the lifespan of aligners, retainers, and dentures while maintaining their function and transparency.



2 Soniprox+ CleanFoam

- ✓ Cleaning foam always at hand, easy to use.
- ✓ Eliminates bacteria and viruses in just a few seconds.
- ✓ Gentle and safe for all materials.



3 Soniprox+ Tab

- ✓ Tablets for deep and regular cleaning.
- ✓ Effectively remove discoloration and deposits.
- ✓ Neutralize unpleasant odors, restoring freshness.



4 Soniprox+ TravelSet

- ✓ An essential hygiene kit for work, school, and travel.
- ✓ Includes a foldable toothbrush and interdental brushes.
- ✓ An elegant, compact case for convenient transport.